

Product Datasheet

Bovine ANGPT1 ELISA Kit (Ready to Use) (orb549160)

Description	Bovine Angiopoietin 1 (ANGPT1) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for co...
Reactivity	Bovine
Range	0.625-40ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	O18920
Sensitivity	0.195ng/mL
Entrez	282140
Expiration Date	Please enquire.