

Product Datasheet

Bovine Angiotensinogen ELISA Kit (Ready to Use) (orb549162)

Description	Bovine Angiotensin II (Angiotensinogen) Sandwich ELISA KitReady to Use: Pre-diluted detection reagen...
Reactivity	Bovine
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at - 20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Sensitivity	8.46pg/mL
Expiration Date	Please enquire.