

Product Datasheet

Bovine AMH ELISA Kit (Ready to Use) (orb549164)

Description	Bovine Anti-Mullerian Hormone (AMH) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents f...
Reactivity	Bovine
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P03972
Sensitivity	0.117ng/mL
Entrez	280718
Expiration Date	Please enquire.