

Product Datasheet

Bovine BMPR1A ELISA Kit (Ready to Use) (orb549175)

Description	Bovine Bone Morphogenetic Protein Receptor 1A (BMPR1A) Sandwich ELISA KitReady to Use: Pre-diluted d...
Reactivity	Bovine
Range	1.56-100ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	F1MCE7
Sensitivity	0.53ng/mL
Entrez	404101
Expiration Date	Please enquire.