

Product Datasheet

Bovine CSN3 ELISA Kit (Ready to Use) (orb549180)

Description	Bovine Casein Kappa (CSN3) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conven...
Reactivity	Bovine
Range	0.312-20µg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitymilk.
Sample Types	milk.
Uniprot ID	P02668
Sensitivity	0.101µg/mL
Entrez	281728
Expiration Date	Please enquire.