

Product Datasheet

Bovine CD40L ELISA Kit (Ready to Use) (orb549186)

Description	Bovine Cluster Of Differentiation 40 Ligand (CD40L) Sandwich ELISA KitReady to Use: Pre-diluted dete...
Reactivity	Bovine
Range	0.625-40ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P51749
Sensitivity	0.18ng/mL
Entrez	282387
Expiration Date	Please enquire.