

Product Datasheet

Bovine CKM ELISA Kit (Ready to Use) (orb549189)

Description	Bovine Creatine Kinase, Muscle (CKM) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents ...
Reactivity	Bovine
Range	12.5-800U/L
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q9XSC6
Sensitivity	4.65U/L
Entrez	286822
Expiration Date	Please enquire.