

Product Datasheet

Bovine DNASE1 ELISA Kit (Ready to Use) (orb549194)

Description	Bovine Deoxyribonuclease I (DNASE1) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents f...
Reactivity	Bovine
Range	0.781-50ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P00639
Sensitivity	0.32ng/mL
Entrez	282217
Expiration Date	Please enquire.