

Product Datasheet

Bovine FGF2 ELISA Kit (Ready to Use) (orb549200)

Description	Bovine Fibroblast Growth Factor 2, Basic (FGF2) Sandwich ELISA KitReady to Use: Pre-diluted detectio...
Reactivity	Bovine
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	F1N3T9
Sensitivity	6.5pg/mL
Entrez	281161
Expiration Date	Please enquire.