

Product Datasheet

Bovine Haptoglobin ELISA Kit (Ready to Use) (orb549212)

Description Bovine Haptoglobin (Haptoglobin) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for ...

Reactivity Bovine

Range 156-10000ng/mL

Form/Appearance Ready to Use

Storage All reagents should be stored at -20°C upon receipt

Note For research use only

Application notes Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.


Sample Types serum, plasma and other biological fluids.

Uniprot ID [Q2TBU0](#)

Sensitivity 41.5ng/mL

Entrez [280692](#)

Expiration Date Please enquire.



Standard
Curve of
Bovine
Haptoglobin
ELI...