

Product Datasheet

Bovine IGFBP6 ELISA Kit (Ready to Use) (orb549235)

Description	Bovine Insulin Like Growth Factor Binding Protein 6 (IGFBP6) Sandwich ELISA KitReady to Use: Pre-dil...
Reactivity	Bovine
Range	2.5-160ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q05718
Sensitivity	0.94ng/mL
Entrez	404186
Expiration Date	Please enquire.