

Product Datasheet

Bovine LPO ELISA Kit (Ready to Use) (orb549248)

Description	Bovine Lactoperoxidase (LPO) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conv...
Reactivity	Bovine
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitymilk, tears or other biological fluids.
Sample Types	milk, tears and other biological fluids.
Uniprot ID	P80025
Sensitivity	0.058ng/mL
Entrez	280844
Expiration Date	Please enquire.