

Product Datasheet

Bovine MMP7 ELISA Kit (Ready to Use) (orb549262)

Description	Bovine Matrix Metalloproteinase 7 (MMP7) Sandwich ELISA KitReady to Use: Pre-diluted detection reage...
Reactivity	Bovine
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q148N3
Sensitivity	0.089ng/mL
Entrez	286794
Expiration Date	Please enquire.