

Product Datasheet

Bovine NGF ELISA Kit (Ready to Use) (orb549266)

Description	Bovine Nerve Growth Factor (NGF) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for ...
Reactivity	Bovine
Range	125-8000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P13600
Sensitivity	30pg/mL
Entrez	281350
Expiration Date	Please enquire.