

Product Datasheet

Bovine Neuropeptide Y ELISA Kit (Ready to Use) (orb549267)

Description	Bovine Neuropeptide Y (Neuropeptide Y) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent...
Reactivity	Bovine
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at - 20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q6RUW3
Sensitivity	8.52pg/mL
Entrez	504216
Expiration Date	Please enquire.