

Product Datasheet

Bovine Osteopontin ELISA Kit (Ready to Use) (orb549273)

Description	Bovine Osteopontin (Osteopontin) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for ...
Reactivity	Bovine
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P31096
Sensitivity	0.069ng/mL
Entrez	281499
Expiration Date	Please enquire.