

Product Datasheet

Bovine PTN ELISA Kit (Ready to Use) (orb549282)

Description	Bovine Pleiotrophin (PTN) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conveni...
Reactivity	Bovine
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P21782
Sensitivity	11.4pg/mL
Entrez	280904
Expiration Date	Please enquire.