

Product Datasheet

Bovine TPO ELISA Kit (Ready to Use) (orb549305)

Description	Bovine Thrombopoietin (TPO) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conve...
Reactivity	Bovine
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	C8CJH1
Sensitivity	31pg/mL
Entrez	525013
Expiration Date	Please enquire.