

Product Datasheet

Canine ADP ELISA Kit (Ready to Use) (orb549327)

Description	Canine Adiponectin (ADP) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenie...
Reactivity	Canine
Range	1.56-100ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q76C76
Sensitivity	0.68ng/mL
Entrez	403625
Expiration Date	Please enquire.