

Product Datasheet

Canine Adrenomedullin ELISA Kit (Ready to Use) (orb549329)

Description	Canine Adrenomedullin (Adrenomedullin) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent...
Reactivity	Canine
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at - 20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	077559
Sensitivity	4.74pg/mL
Entrez	403817
Expiration Date	Please enquire.