

Product Datasheet

Canine aFP ELISA Kit (Ready to Use) (orb549330)

Description	Canine Alpha-Fetoprotein (aFP) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for co...
Reactivity	Canine
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q8MJU5
Sensitivity	0.133ng/mL
Entrez	403551
Expiration Date	Please enquire.