

Product Datasheet

Canine CLU ELISA Kit (Ready to Use) (orb549344)

Description	Canine Clusterin (CLU) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenient...
Reactivity	Canine
Range	1.56-100ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P25473
Sensitivity	0.5ng/mL
Entrez	442971
Expiration Date	Please enquire.