

Product Datasheet

Canine factor viii ELISA Kit (Ready to Use) (orb549345)

Description	Canine Coagulation Factor VIII (factor viii) Sandwich ELISA KitReady to Use: Pre-diluted detection r...
Reactivity	Canine
Range	3.12-200ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityplasma or other biological fluids.
Sample Types	plasma and other biological fluids.
Uniprot ID	O18806
Sensitivity	1.15ng/mL
Entrez	403875
Expiration Date	Please enquire.