

Product Datasheet

Canine CX43 ELISA Kit (Ready to Use) (orb549349)

Description	Canine Connexin 43 (CX43) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conveni...
Reactivity	Canine
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	J9P7X9
Sensitivity	0.074ng/mL
Entrez	403418
Expiration Date	Please enquire.