

Product Datasheet

Canine MCP2 ELISA Kit (Ready to Use) (orb549387)

Description	Canine Monocyte Chemotactic Protein 2 (MCP2) Sandwich ELISA KitReady to Use: Pre-diluted detection r...
Reactivity	Canine
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q68AY9
Sensitivity	9.7pg/mL
Entrez	448792
Expiration Date	Please enquire.