

Product Datasheet

Canine TPO ELISA Kit (Ready to Use) (orb549400)

Description	Canine Thrombopoietin (TPO) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conve...
Reactivity	Canine
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P42705
Sensitivity	6.9pg/mL
Entrez	100855822
Expiration Date	Please enquire.