

Product Datasheet

Gallus BMP4 ELISA Kit (Ready to Use) (orb549414)

Description	Chicken Bone Morphogenetic Protein 4 (BMP4) Sandwich ELISA KitReady to Use: Pre-diluted detection re...
Reactivity	Gallus
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	E1C0T2
Sensitivity	12.6pg/mL
Entrez	396165
Expiration Date	Please enquire.