

## Product Datasheet

### Gallus CNTF ELISA Kit (Ready to Use) (orb549420)

|                          |                                                                                                                                                           |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | Chicken Ciliary Neurotrophic Factor (CNTF) Sandwich ELISA KitReady to Use: Pre-diluted detection rea...                                                   |
| <b>Reactivity</b>        | Gallus                                                                                                                                                    |
| <b>Range</b>             | 62.5-4000pg/mL                                                                                                                                            |
| <b>Form/Appearance</b>   | Ready to Use                                                                                                                                              |
| <b>Storage</b>           | All reagents should be stored at -20°C upon receipt                                                                                                       |
| <b>Note</b>              | For research use only                                                                                                                                     |
| <b>Application notes</b> | Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.                                                   |
| <b>Uniprot ID</b>        | <b>Q02011</b>                                                                                                                                             |
| <b>Sensitivity</b>       | 29.4pg/mL                                                                                                                                                 |
| <b>Entrez</b>            | <b>396488</b>                                                                                                                                             |
| <b>Expiration Date</b>   | Please enquire.                                                                                                                                           |