

Product Datasheet

Gallus FGF9 ELISA Kit (Ready to Use) (orb549428)

Description	Chicken Fibroblast Growth Factor 9 (FGF9) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent...
Reactivity	Gallus
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q7ZZN4
Sensitivity	35pg/mL
Entrez	378917
Expiration Date	Please enquire.