

Product Datasheet

Gallus IL12B ELISA Kit (Ready to Use) (orb549442)

Description	Chicken Interleukin 12B (IL12B) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for c...
Reactivity	Gallus
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Sensitivity	7.2pg/mL
Expiration Date	Please enquire.