

Product Datasheet

Gallus IL16 ELISA Kit (Ready to Use) (orb549444)

Description	Chicken Interleukin 16 (IL16) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for con...
Reactivity	Gallus
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q70XC4
Sensitivity	5.7pg/mL
Entrez	374270
Expiration Date	Please enquire.