

Product Datasheet

Gallus LEP ELISA Kit (Ready to Use) (orb549450)

Description	Chicken Leptin (LEP) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenient a...
Reactivity	Gallus
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	O42164
Sensitivity	14.8pg/mL
Entrez	373955
Expiration Date	Please enquire.