

Product Datasheet

Equine DEFb2 ELISA Kit (Ready to Use) (orb549469)

Description	Equine Defensin Beta 2 (DEFb2) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for co...
Reactivity	Equine
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q0W9Q1
Sensitivity	28pg/mL
Entrez	100629494
Expiration Date	Please enquire.