

Product Datasheet

Equine ECF ELISA Kit (Ready to Use) (orb549470)

Description	Equine Eosinophil Chemotactic Factor (ECF) Sandwich ELISA KitReady to Use: Pre-diluted detection rea...
Reactivity	Equine
Range	62.5-4000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9TTQ4
Sensitivity	26.1pg/mL
Entrez	100033949
Expiration Date	Please enquire.