

Product Datasheet

Equine MMP1 ELISA Kit (Ready to Use) (orb549490)

Description	Equine Matrix Metalloproteinase 1 (MMP1) Sandwich ELISA KitReady to Use: Pre-diluted detection reage...
Reactivity	Equine
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9XSZ5
Sensitivity	0.040ng/mL
Entrez	100033896
Expiration Date	Please enquire.