

Product Datasheet

Equine TGF beta 1 ELISA Kit (Ready to Use) (orb549501)

Description	Equine Transforming Growth Factor Beta 1 (TGF beta 1) Sandwich ELISA KitReady to Use: Pre-diluted de...
Reactivity	Equine
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, platelet-poor plasma, tissue homogenates, cell culture supernatants and other biological fluids.
Uniprot ID	O19011
Sensitivity	6.3pg/mL
Entrez	100033900
Expiration Date	Please enquire.