

Product Datasheet

Equine TGFb3 ELISA Kit (Ready to Use) (orb549502)

Description	Equine Transforming Growth Factor Beta 3 (TGFb3) Sandwich ELISA KitReady to Use: Pre-diluted detecti...
Reactivity	Equine
Range	125-8000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	F7DR45
Sensitivity	31pg/mL
Entrez	100051765
Expiration Date	Please enquire.