

Product Datasheet

Equine TPS ELISA Kit (Ready to Use) (orb549503)

Description	Equine Tryptase (TPS) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenient ...
Reactivity	Equine
Range	0.625-40ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q7YS62
Sensitivity	0.229ng/mL
Entrez	100033966
Expiration Date	Please enquire.