

Product Datasheet

Other HVD3 ELISA Kit (Ready to Use) (orb549515)

Description	General 25-Hydroxyvitamin D3 (HVD3) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents f...
Reactivity	All
Range	3.12-200ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Sensitivity	0.83ng/mL
Expiration Date	Please enquire.