

Product Datasheet

Other ALD ELISA Kit (orb549524)

Description	General Aldosterone (ALD) Sandwich ELISA Kit
Reactivity	All
Range	31.2-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Sensitivity	11.2pg/mL
Expiration Date	Please enquire.