

Product Datasheet

Other HA ELISA Kit (Ready to Use) (orb549566)

Description	General Hyaluronic Acid (HA) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conv...
Reactivity	All
Range	6.25-400ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Sensitivity	2.17ng/mL
Expiration Date	Please enquire.