

Product Datasheet

Other NAADP ELISA Kit (Ready to Use) (orb549573)

Description	General Nicotinic Acid Adenine Dinucleotide Phosphate (NAADP) Sandwich ELISA KitReady to Use: Pre- di...
Reactivity	All
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Sensitivity	0.055ng/mL
Expiration Date	Please enquire.