

Product Datasheet

Other VB6 ELISA Kit (Ready to Use) (orb549600)

Description	General Vitamin B6 (VB6) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenie...
Reactivity	All
Range	3.12-200ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Sensitivity	1.17ng/mL
Expiration Date	Please enquire.