

Product Datasheet

Goat BMP4 ELISA Kit (Ready to Use) (orb549606)

Description	Goat Bone Morphogenetic Protein 4 (BMP4) Sandwich ELISA KitReady to Use: Pre-diluted detection reage...
Reactivity	Goat
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	A6N998
Sensitivity	0.148ng/mL
Entrez	100860789
Expiration Date	Please enquire.