

Product Datasheet

Goat GHRH ELISA Kit (Ready to Use) (orb549608)

Description	Goat Growth Hormone Releasing Hormone (GHRH) Sandwich ELISA KitReady to Use: Pre-diluted detection r...
Reactivity	Goat
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P63293
Sensitivity	6.4pg/mL
Entrez	102181610
Expiration Date	Please enquire.