

Product Datasheet

Guinea pig C3 ELISA Kit (Ready to Use) (orb549620)

Description	Guinea pig Complement Component 3 (C3) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent...
Reactivity	Guinea pig
Range	1.56-100ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Sensitivity	0.64ng/mL
Expiration Date	Please enquire.