

Product Datasheet

Human OXCT1 ELISA Kit (Ready to Use) (orb549640)

Description	Human 3-Oxoacid Coenzyme A Transferase 1 (OXCT1) Sandwich ELISA KitReady to Use: Pre-diluted detecti...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P55809
Sensitivity	0.124ng/mL
Entrez	5019
Expiration Date	Please enquire.