

Product Datasheet

Human ADAM19 ELISA Kit (Ready to Use) (orb549656)

Description	Human A Disintegrin And Metalloprotease 19 (ADAM19) Sandwich ELISA KitReady to Use: Pre-diluted dete...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9H013
Sensitivity	0.108ng/mL
Entrez	8728
Expiration Date	Please enquire.