

Product Datasheet

Human ABI1 ELISA Kit (Ready to Use) (orb549668)

Description	Human Abl Interactor 1 (ABI1) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for con...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q8IZP0
Sensitivity	0.055ng/mL
Entrez	10006
Expiration Date	Please enquire.