

Product Datasheet

Human ACTb ELISA Kit (Ready to Use) (orb549687)

Description	Human Actin Beta (ACTb) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenien...
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P60709
Sensitivity	0.63ng/mL
Entrez	60
Expiration Date	Please enquire.